boso medistar+

**BOSCH
+SOHN**

boso

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2

**User
Instructions**

42

**Mode
d'emploi**

82

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122

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162

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Table of contents

Product features	44
Contents of package	45
Symbols on the blood pressure monitor	46
Quick guide	49
Introduction	50
Purpose	52
Blood pressure values	53
Systolic and diastolic blood pressure value	53
WHO blood pressure cut-off values	53
Starting up	55
Inserting batteries	55
Battery status	56
Changing batteries	56
General instructions for self-measurement	57
Preparing to measure your blood pressure	59

Note on cuff sizes.....	59
How to wear the unit.....	60
Measuring your blood pressure	63
Interrupting a measurement	64
Measurement display	65
Completing the blood pressure record card / Retrieving the average value and classification scale.....	67
Displaying individual stored values	70
Error messages	71
Cleaning and Disinfection.....	73
Warranty Conditions	74
Customer service / Disposal.....	75
Technical data.....	76
Calibration checks - Testing instructions (for trained specialist staff only)	79

44 Product features

EN

Battery compartment

WHO
classification

LCD display

M button

START button

Contents of package

1 Blood pressure monitor
boso-
medistar+

2 LR 03 (AAA)

1 Case


1 User instructions

1 Blood pressure record card

Symbols on the blood pressure monitor

Symbol	Function/meaning
	START/STOP button
 	M button (Memory button)
	Serial number
	Year of manufacture
	Type BF: The device, cuff and tubing have been designed to protect the patient against electric shock.

Symbols on the blood pressure monitor

Symbol	Function/meaning
	Reading in progress. Flashes as soon as a heartbeat was detected.
	Cardiac rhythm disorder or movement during reading.
M	Stored value
	Battery fully charged.
	Replace battery as soon as possible.
	Flashing! No further reading possible. Replace battery.
	Protect against liquids

Symbols on the blood pressure monitor

Symbol	Function/meaning
<i>SYS.</i>	Systolic blood pressure in mmHg
<i>DIA.</i>	Diastolic blood pressure in mmHg
<i>PUL.</i>	Pulse per minute
 CE 0124	Device complies with the European Medical Devices Directive.
	Device must not be discarded with household waste.
	Manufacturer
	Read instructions for use

Quick guide

➔ Insert batteries, ensure correct positioning (page 55)

➔ Attach unit (page 60)

 When measuring, the wrist must be held at the same height as the heart!

➔ Start (and stop) a measurement by pressing the 'START' button (page 63)

Measured systolic, diastolic and pulse values are displayed after the measurement is completed (page 65)

➔ Retrieving values from memory (page 67): with the instrument switched off, press the 'M' button. The

first value displayed is the average value of all stored measurements. The total number of stored measurements is displayed at the same time, e.g. 'A90'. After repeated pressing the M button, each individual measurement is displayed.

 This quick guide does not replace the detailed information on the use and safety of the blood pressure monitor featured in the user instructions.

Therefore, please ensure reading the user instructions!

Dear Customer,

We are delighted that you have decided to purchase a bosco blood pressure monitor. The bosco brand is synonymous with optimal quality and precision and is also number 1 among professionals: 96% of all German doctors currently use bosco blood pressure monitors (GfK 1/2010).

The experience gained from decades of use by professionals has also been incorporated into all the patient monitors for home measurement.

This instrument has also passed our strict quality control procedures and is

your reliable partner for monitoring your blood pressure values.

 Please read these user instructions before using the instrument for the first time, as correct blood pressure readings can only be obtained if the instrument is operated correctly.

In these user instructions, '➔' indicates an action by the user.

To get help concerning usage or maintenance please contact your dealer or the manufacturer (contact details can be found on the back cover of this manual).

The bosomedistar+ blood pressure monitor is suitable for patients of all ages with a wrist circumference of between 13,5 and 21,5 cm.

The use of accessories not referred to in these instructions for use can undermine safety.

Wireless communication devices, such as home networking devices, mobile phones, cordless phones and their base stations, walkie-talkies can affect this blood pressure monitor. Therefore, a minimum distance of 11 feet should be kept from such devices.

An unexpected operating condition or an incidence which has worsened the health condition or could have worsened the health condition should be reported to the manufacturer immediately.

If the monitor is used in medical practice (in accordance with the legislation governing the operation of medical devices), it should be calibrated at regular intervals (see page 79).

If the instrument is to be sold, ensure that these user instructions are enclosed.

Purpose

Non-invasive recording of systolic and diastolic blood pressure and pulse frequency in humans.

Conditions or diseases to be verified, monitored, treated or diagnosed:

- hypertension
- arrhythmia

Not suitable for use on newborn babies and infants.

The device is suitable for self-measurement at home. The patient is regarded as the operator.

Blood pressure values

To determine the blood pressure, **two values** need to be measured:

- **Systolic (upper) blood pressure:**
This is produced when the heart muscle contracts and blood is forced into the blood vessels.

- **Diastolic (lower) blood pressure:**
This is present when the heart muscle expands and fills with blood again.

Blood pressure measurement values are expressed in mmHg (mm of mercury).

The World Health Organization (WHO, 1999) has defined the following cut-off values for assessing blood pressure:

	systolic	diastolic
High blood pressure	over 140 mmHg	over 90 mmHg
Borderline normal blood pressure	130 to 139 mmHg	85 to 89 mmHg
Normal blood pressure	120 to 129 mmHg	80 to 84 mmHg
Optimum blood pressure	up to 119 mmHg	up to 79 mmHg

Blood pressure values

Whether medical treatment is necessary depends not only on the blood pressure, but also on the patient's risk profile. Please contact your GP if one of the values (systole, diastole) consistently exceeds the cut-off point for high blood pressure.

The logo consists of the letters 'EN' in a white, sans-serif font, enclosed within a white circle. This circle is centered on a dark grey rectangular background.

Self-measurement of blood pressure often produces somewhat lower values at home than at the doctor's.

The German Hypertension League has therefore defined lower cut-off values for home measurements:

Home measurement: 135/85 mmHg

Measurement in the
doctor's surgery: 140/90 mmHg

Starting up

The device must be installed and put into operation in accordance with the instructions in this manual.


Illustration 1

- ➔ Use only high quality leak-proof batteries of the correct specifications (see 'Technical data' page 76).

⚠ **Never** mix old and new batteries or different makes.

⚠ If the batteries are not inserted correctly, the monitor will not work and this may cause the batteries to heat up and leak, damaging the monitor.

➔ **Inserting the batteries**

The battery compartment is located on the left side of the monitor. Insert batteries in the battery compartment as shown in illustration 1.

- ➔ Remove the batteries if the monitor is not to be used for any length of time.

⚠ Batteries that are short-circuited can become hot and cause burning.

Starting up

The monitor has a battery status display.

Batteries fully charged.

Replace batteries soon.

Flashing! No further measurements possible.
Replace batteries.

EN

⚠ Change of battery deletes memory!

The removal of batteries will erase the memory store immediately.

⚠ Do not touch the batteries and the patient at the same time.

⚠ In the case of battery short circuit, the temperature of the cuff may get hot up to 50°C and the temperature of the complete unit may become hot. If you feel high temperature, remove the cuff immediately from the wrist.

⚠ Care for the environment

Do not dispose of used batteries and charger units in domestic waste. You can take these to collection sites for used batteries or special waste. Contact your local authority for information.

1. Fluctuations in blood pressure are quite normal. Even when measurements are repeated, there can be marked differences. Single or irregular measurements do not provide reliable information about the actual blood pressure. A reliable evaluation is only possible if regular measurements are taken under comparable conditions and the measurement values are entered in the blood pressure record card.
2.  Self-testing does not mean self-treatment. Do not on any account change the dose of medicine prescribed by your doctor of your own accord.
3. Irregular heartbeats can affect the accuracy of the monitor or result in incorrect measurements (see page 72).
4. Incorrect measurements can also occur in patients with a weak pulse (possibly in patients wearing a pacemaker, for example). The blood pressure monitor itself has no effect on the pacemaker.
5. The device has not been validated for use on pregnant women.

6. Always take the measurement in a calm and relaxed state.

You are recommended to check your blood pressure twice a day – in the morning on getting up and in the evening once you have relaxed after work.

A small grey square icon containing the letters 'EN' in white, enclosed in a white circle.

7. Always measure your blood pressure from the wrist with the higher blood pressure values (where not required otherwise (see page 61)).

To do this, measure the blood pressure first in both wrists and thereafter always from the wrist with the higher blood pressure.

Note on cuff sizes:

The unit must only be used in combination with the CW 61 cuff (wrist circumference 13.5 to 21.5 cm).

Observe the following instructions when measuring your blood pressure:

Avoid smoking or drinking coffee for one hour before measuring your blood pressure.

Sit comfortably for the blood pressure measurement. Support your back and arms. Do not cross your legs. Place your feet flat on the floor. Relax for 5 minutes before taking a measurement. Do not move during the measurement.

Preparing to measure your blood pressure

How to wear the unit

The measurement should be taken from the bare wrist.

 Please make sure the display is positioned at your inner wrist (see illustration 2).


Illustration 2

➔ Close the cuff with the velcro fastener. The cuff should fit snugly around the wrist. When relaxed, place your elbow on your chest with the wrist in the same height as your heart (see illustration 3).


Illustration 3

 Do not attach the cuff over wounds as this can cause further injuries.

 Ensure that the cuff is not attached to an arm where arteries or veins are undergoing or have undergone medical treatment (e.g. shunt).

 For women who have had a breast amputated, do not attach the cuff to the wrist on the amputated side of the body.

 During the measurement, medical apparatus that are being used on the same arm at the same time, may malfunction.

 It can occur that there are differences between measurements at the wrist and the upper arm.

Therefore, it is recommended on occasions to take comparative measurements at the upper arm.

 In cases of poor blood circulation or severe arteriosclerosis, it is not appropriate to take measurements at the wrist. In these instances the blood pressure should be taken at the upper arm.

 **When measuring, the wrist must be held at the same height as the heart!**

Preparing to measure your blood pressure

 The blood pressure measurement should not prevent the flow of blood for an unnecessary length of time (> 2 minutes). If the monitor fails to function correctly, remove the unit from the wrist.

EN

 Excessively frequent measurements can cause injury by restricting blood flow.

 The device may not be used by unsupervised children.

 Do not use the device near infants. This can lead to accidents or damage.

 Do not start the device without putting on the cuff.

 There are small parts that may cause a choking hazard if swallowed by mistake by infants.

 If the ambient temperature changes significantly, leave the device about an hour to adapt to the new ambient temperature before using it.

 The performance of the device can be affected by excessive temperature, humidity or altitude.

Measuring your blood pressure

Your monitor has a memory with a capacity of 90 measurements.

- ➔ Start the measurement with the 'START' button

- ⚠ **Now hold the arm absolutely still and do not talk.**

All the items on the LCD display appear briefly to indicate that the monitor is now ready for use (see illustration 4).

Illustration 4

Measuring your blood pressure

The pump starts to inflate the cuff.

The display indicates the pressure of the cuff as it rises.

The unit inflates the cuff to a pressure of about 200 mmHg. The display indicates the pressure of the cuff as it rises. In case this pressure is not sufficient, the unit inflates the cuff by a further 50 mmHg.

Once inflated, the pump is switched off and the air in the cuff escapes automatically.

If not already the case, you should now keep your arm completely still and not speak.

As soon as the unit registers the first beats, the '♥' symbol blinks in time with the pulse rhythm.

Interrupting a measurement:

A measurement can be interrupted at any stage by pressing the 'START' button; the cuff deflates automatically.

Measurement display

After the end of the measurement, the inbuilt valve opens automatically and the cuff deflates rapidly.

The measured blood pressure values (systole, diastole, pulse) are displayed (see illustration 5).

III. 5

The measurement is stored automatically. If the memory store is full (90 measurements), the oldest measurement is deleted and the current one stored.

Guest mode:

If a measurement value is not to be kept, press „M” while the measured values are still on display.

Measurement display

Invalid measurements with an 'Err' message (see page 71) are not stored.

If the symbol '❤️' appears after the measurement, it is recommended that the measurement is repeated, holding the arm completely still.

If the '❤️' symbol appears again even with the arm held still, this may be an indication of an irregular heartbeat. Please discuss this with your doctor at your next visit.

Blood pressure is a dynamic parameter and can be affected by the patient's position, e.g. sitting, standing, lying, moving, before or during the measurement as well as by his physical condition (stress, disease, etc.).

- ➔ If the measurement result is obviously wrong, repeat the measurement.

The unit switches itself off automatically after about 1 min. To repeat a measurement, press the 'START' button again.

Allow at least two minutes between measurements.

- ➔ If no further measurement is required, remove the unit from the wrist.

Retrieving the average value and classification scale

- ➔ Enter each individual value in your record card.
- ➔ Enter the average value in the field specifically reserved for this in your record card.
- ➔ Recall the individual average value after 30 measurements by pressing the 'M' button while the unit is switched off. The average value of all stored values (systolic and diastolic) appears. The display <pulse> shows the number of stored measurements (see illustration 6).

If there is no reading in the memory, "--" will be displayed for systole, diastole and pulse.

Retrieving the average value and WHO classification

A dark rectangle appears on the left of the screen (see illustration 6), indicating the level of the average blood pressure value.

EN


Illustration 6

Retrieving the average value and WHO classification

Use the accompanying evaluation scale (illustration 7) to identify the range within which your blood pressure falls under the WHO classification (see page 53).

EN

WHO classification

-  High blood pressure (hypertension)
-  Borderline normal blood pressure
-  Normal blood pressure
-  Optimum blood pressure

Illustration 7

Displaying individual stored values

- ➔ In order to bring back the last measurement stored, press the M button once more. The number of the measurements appears in the display (see illustration 8). After two seconds the corresponding measurement value (systole, diastole and pulse) appears automatically (see illustration 9).

Illustration 8

Illustration 9

Systole
mmHg
Diastole
mmHg
Pulse/min.



- ➔ By repeated pressing of the M button, all stored measurements will be brought up in succession and appear in the display as described above.

If no button is pressed, the unit switches off automatically after 5 seconds.

- ➔ To start a new measurement, simply press „START“.
- ➔ **How to erase memory:**
Push M button during memory display for approx. 5 seconds until „M“ starts blinking in the display.

If any problems occur during the measurement, the display will show an error message instead of a blood pressure reading.

Meaning of error messages:

Cause of Err message	Troubleshooting
 Err 1	Keep arm still Repeat measurement.
 Err 2	Check the positioning of the cuff. Repeat measurement.
 Err CUF	Cuff is possibly not tight enough. Repeat measurement.

Cause of Err message	Troubleshooting
 <p>The pulse value could not be determined correctly.</p>	<p>Check the positioning of the cuff. Repeat measurement.</p>
 <p>Irregular puls or incorrect taking of blood pressure (i.e. movement during measurement)</p>	<p>Repeat measurement holding the arm absolutely still. If the “♥” symbol is displayed again, it can be that an irregular heart beat is present. Consult your GP.</p>

Cleaning Monitor:

Use a soft, dry cloth to clean your monitor.

Cleaning Cuff:

Small stains on the cuff can be removed carefully with a proprietary washing-up liquid.

Disinfection Cuff:

For disinfectant wipes (at least 5 minutes exposure time) of the device and the cuff, we recommend the disinfectant mikro-zid sensitive liquid (Schülke & Mayr). To disinfect the cuff, we recommend spray disinfection. In particular, if the device is used on multiple users is to pay attention to regular cleaning and disinfections of the cuff.

We give 3 years warranty from the date of purchase. The purchase date has to be proven by the invoice. Within the warranty period defects are eliminated free of charge. After repairs the warranty period is not extended on the whole unit, but only to the replaced components.

EN

Excluded from the warranty are parts subject to normal wear and tear (e.g. cuff), transport damages and any damage caused by improper handling (e.g. non-compliance with the instructions for use). Damages due to disassembly by unauthorized persons are also excluded from warranty.

No claims for damages against us are substantiated by the warranty.

In the case of justified warranty claims the device has to be sent along with the original invoice to:

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Warranty and repair work must be undertaken by trained and authorised personnel. Do not modify this equipment without authorization of the manufacturer.

Send the monitor, carefully packaged and with sufficient postage, to your authorised dealer or directly to:

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Do not dispose of monitors and batteries in the domestic waste.

At the end of its lifespan, the monitor must be taken to a collection site for obsolete electronic items.

Please note that batteries and rechargeable batteries must be disposed of separately (local authority collection site).

Expected operational lifetime of the equipment: 10 years.

Expected operational lifetime of the cuff: 10.000 measurement cycles.

Technical data

Measurement principle:	oscillometric
Measurement range:	20 to 280 mmHg, 40 to 180 pulses/min.
Cuff pressure:	0 to 299 mmHg
 Memory store:	90 measurements
Display:	LCD
Operating conditions:	environmental temperature +10°C to +40°C relative humidity 15 to 85% Air pressure 800 hPa to 1060 hPa
Transport/ Storage conditions:	environmental temperature -20°C to +60°C relative humidity 10 to 95% Air pressure 700 hPa to 1060 hPa

Power supply:	DC 3 V (2 x 1.5 V Micro IEC LR 03 alkaline manganese batteries)
Typical battery life:	250 measurement cycles (depending on inflation pressure and frequency of use).
Battery check:	symbol display in LCD window
Weight:	100 g without batteries
Dimensions (WxHxD):	65 mm x 70 mm x 70 mm
Classification:	Medical device with internal energy source, continuous operation mode
Protection against solid objects:	IP20 IP classification is the degree of protection provided by enclosures in accordance with IEC 60529. This device is protected against solid foreign objects of 12mm diameter and greater such as fingers. This device is not protected against water.

Technical data

**Clinical test
(DIN 58130):**

accuracy complies with the requirements
of EN 1060 Part 3

**Maximum deviation of cuff
pressure measurement:**

± 3 mmHg or 2% of the reading
(whichever is the higher)

**Maximum deviation of
pulse rate display:**

± 5 %

Applicable standard:

IEC 80601-2-30 : "Particular requirements for
basic safety and essential performance of
automated non-invasive sphygmomanometers"

(for trained specialist staff only)

Calibration checks – every 2 years at the latest – must be performed by the following institutions or persons:

- Manufacturer
- Trade measurements authorities
- Persons who fulfil the requirements of the legislation governing the operation of medical devices.

A) Functional testing

Functional testing of the unit can only be performed on human subjects or using a suitable simulator.

B) Leak testing of the pressure circuit and error testing of the display

Please note:

- a) Prior to testing, separate the cuff from the unit (after having first removed the screws in the battery compartment). The air hose connecting plug (special plug) can now be inserted into the unit.

(for trained specialist staff only)

b) To perform leak testing of the pressure circuit, a tee connector must be used to connect the rolled-up cuff and a pump bulb to the pressure circuit.

c) If the pressure in measuring mode exceeds 320 mmHg, the quick-action air release valve triggers and opens the pressure circuit. If the pressure in testing mode exceeds 320 mmHg, the display blinks.

Testing

1. Remove batteries.
2. Hold START-button pressed down and insert batteries.
3. Release START-button.

The unit is now ready for testing. The displays of SYS and DIA now show the actual pressure.

(for trained specialist staff only)

4. Testing for divergence of Pressure indication and Density of Pressure circuit (setting time for cuff at least 60 seconds) can now be carried out in the usual manner.

5. Switch off unit and reconnect the cuff to the unit.

C) Security

For security purposes the two parts of the housing are joined together with a security label.

 The device must not undergo service or maintenance while it is in use.

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BOSCH + SOHN GmbH u. Co. KG